

SPATA2L 抗原（重组蛋白）

中文名称： SPATA2L 抗原（重组蛋白）

英文名称： SPATA2L Antigen (Recombinant Protein)

别名： tamo, C16orf76

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 153-424 amino acids of human SPATA2L

技术规格

Full name:	spermatogenesis associated 2-like
Synonyms:	tamo, C16orf76
Swissprot:	Q8IUW3
Gene Accession:	BC039608
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene belongs to the SPATA2 family. High levels expression in pancreas.