

## 兔抗 ERP29 多克隆抗体

中文名称：兔抗 ERP29 多克隆抗体

英文名称：Anti-ERP29 rabbit polyclonal antibody

别名：ERp28; ERp31; PDIA9; PDI-DB; C12orf8

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ERP29

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a reticuloplasmin, a protein which resides in the lumen of the endoplasmic reticulum (ER). The protein shows sequence similarity to the protein disulfide isomerase family. However, it lacks the thioredoxin motif characteristic of this family, suggesting that this protein does not function as a disulfide isomerase. The protein dimerizes and is thought to play a role in the processing of secretory proteins within the ER. Alternative splicing results in multiple transcript variants encoding different isoforms
<b>Applications:</b>	ELISA, WB, IHC

<b>Name of antibody:</b>	ERP29
<b>Immunogen:</b>	Synthetic peptide of human ERP29
<b>Full name:</b>	endoplasmic reticulum protein 29
<b>Synonyms :</b>	ERp28; ERp31; PDIA9; PDI-DB; C12orf8
<b>SwissProt:</b>	P30040
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human colon cancer and Human liver cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	29 kDa
<b>WB Positive control:</b>	Hela cells and human fetal liver tissue, A549 and 293T c ells
<b>WB Recommended dilution:</b>	500-2000



