

兔抗 ESR1(Ab-537) 多克隆抗体

中文名称: 兔抗 ESR1(Ab-537) 多克隆抗体

英文名称: Anti-ESR1(Ab-537) rabbit polyclonal antibody

别名: ER; ESR; Era; ESRA; ESTRR; NR3A1

抗原: ESR1(Ab-537)

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

反应种属: Human Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	ER; ESR; Era; ESRA; ESTRR; NR3A1
Applications:	WB, IHC
Name of antibody:	ESR1(Ab-537)
Immunogen:	Synthesized non-phosphopeptide derived from human Estrogen Receptor- α around the phosphorylation site of tyrosine 537 (P-L-Y(p)-D-L).
Full name:	estrogen receptor 1
Synonyms :	ER; ESR; Era; ESRA; ESTRR; NR3A1
SwissProt:	P03372
IHC positive control:	Human lung carcinoma tissue
IHC Recommend dilution:	50-100

WB Predicted band size:	66 kDa
WB Positive control:	293 cells lysate
WB Recommended dilution:	500-3000



