

SYT3 抗原(重组蛋白)

中文名称: SYT3 抗原(重组蛋白)

英文名称: SYT3 Antigen (Recombinant Protein)

别 名: SytIII

储 存: 冷冻(-20℃)

相关类别: 抗原

概述

Fusion protein corresponding to a region derived from 76-200 amino acids of human SYT3

技术规格

Full name:	synaptotagmin III
Synonyms:	SytIII
Swissprot:	Q9BQG1
Gene Accession:	BC031067
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Synaptotagmins are a large gene family of synaptic vesicle type III i ntegral membrane proteins that function as regulators of both exoc ytosis and endocytosis and are involved in neurotransmitter secretion from small secretory vesicles. Calcium binds to Synaptotagmin I which triggers neurotransmitter release at the synapse. Synaptotagmin II is phosphorylated by WNK1 in a process that regulates calcium-dependent interactions. Synaptotagmin III is involved in calcium-dependent



ndent exocytosis of secretory vesicles in endocrine cells and neuron s. Synaptotagmin IV is expressed in neuronal tissues, and has the hi ghest mRNA levels in the hippocampus. The proximity of the Synap totagmin IV gene to markers of several psychiatric disorders sugges t an involvement of synaptotagmin IV in human disease. Synaptotag min V is a dense-core vesicle-specific protein that regulates a specific type of calcium-regulated secretion. Synaptotagmin VI interacts wi th adaptor protein-2 in a calcium-independent manner. Synaptotag min VII is widely expressed in non-neuronal tissues.