

兔抗 ESRRB 多克隆抗体

中文名称：兔抗 ESRRB 多克隆抗体

英文名称：Anti-ESRRB rabbit polyclonal antibody

别名：ERR2; ERRb; ESRL2; NR3B2; DFNB35

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ESRRB

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a protein with similarity to the estrogen receptor. Its function is unknown; however, a similar protein in mouse plays an essential role in placental development.
Applications:	ELISA, IHC
Name of antibody:	ESRRB
Immunogen:	Synthetic peptide of human ESRRB
Full name:	estrogen-related receptor beta
Synonyms :	ERR2; ERRb; ESRL2; NR3B2; DFNB35
SwissProt:	O95718
ELISA Recommended dilution:	2000-5000

IHC positive control:	Human thyroid cancer and Human ovarian cancer
IHC Recommend dilution:	50-200

