

UTP11 抗原（重组蛋白）

中文名称： UTP11 抗原（重组蛋白）

英文名称： UTP11 Antigen (Recombinant Protein)

别名： CGI94; CGI-94; UTP11L

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

Full name:	UTP11 small subunit processome component
Synonyms:	CGI94; CGI-94; UTP11L
Swissprot:	Q9Y3A2
Gene Accession:	BC005182
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Involved in nucleolar processing of pre-18S ribosomal RNA.