

UTY 抗原(重组蛋白)

中文名称: UTY 抗原(重组蛋白)

英文名称: UTY Antigen (Recombinant Protein)

别名: UTY1; KDM6AL

储 存: 冷冻(-20℃)

相关类别: 抗原

概述

Full length fusion protein

技术规格

Full name:	ubiquitously transcribed tetratricopeptide repeat containing , Y-linked
Synonyms:	UTY1; KDM6AL
Swissprot:	O14607
Gene Accession:	BC071744
Purity:	>85%, as determined by Coomassie blue stained SDS-PAG
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a protein containing tetratricopeptide r epeats which are thought to be involved in protein-protein in interactions. The encoded protein is also a minor histoc ompatibility antigen which may induce graft rejection of male stem cell grafts. A large number of alternatively splic



ed transcripts have been observed for this gene, but the f
ull length nature of some of these variants has not been
determined.