

ZEB2 抗原(重组蛋白)

中文名称: ZEB2 抗原(重组蛋白)

英文名称: ZEB2 Antigen (Recombinant Protein)

别 名: SIP1; SIP-1; ZFHX1B; HSPC082; SMADIP1

储 存: 冷冻(-20℃)

相关类别: 抗原

概述

Fusion protein corresponding to a region derived from 1084-1214 amino acids of human ZEB2

技术规格

Full name:	zinc finger E-box binding homeobox 2
Synonyms:	SIP1; SIP-1; ZFHX1B; HSPC082; SMADIP1
Swissprot:	O60315
Gene Accession:	BC127102
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a member of the Zfh1 family of 2-handed zinc finger/homeodomain proteins. It is located in the nucleus and functions as a DNA-binding transcriptional repressor th at interacts with activated SMADs. Mutations in this gene are assoc iated with Hirschsprung disease/Mowat-Wilson syndrome. Alternatively spliced transcript variants have been found for this gene.