

ABHD1 抗原（重组蛋白）

中文名称：ABHD1 抗原（重组蛋白）

英文名称：ABHD1 Antigen (Recombinant Protein)

别名：abhydrolase domain containing 1; LABH1

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to C terminal 300 amino acids of human ABHD1

技术规格：

Full name:	abhydrolase domain containing 1
Synonyms:	LABH1
Swissprot:	Q96SE0
Gene Accession:	BC039576
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene is a member of the AB hydrolase superfamily and encodes a protein with an alpha/beta hydrolase fold. This domain is common to a number of hydrolytic enzymes of widely differing phylogenetic origins and catalytic functions.