

ACTL8 抗原（重组蛋白）

中文名称：ACTL8 抗原（重组蛋白）

英文名称：ACTL8 Antigen (Recombinant Protein)

别名：Actin-like 8; CT57

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to C terminal 250 amino acids of human ACTL8

技术规格：

Full name:	Actin-like 8
Synonyms:	CT57
Swissprot:	Q9H568
Gene Accession:	BC028909
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This protein belongs to the actin family. Actin-like 8 is a protein in humans that is encoded by the ACTL8 gene. The high level expression in testis and pancreas and weak expression in placenta. The subcellular location of this protein in cytoplasm, cytoskeleton.