

ACTR2 抗原(重组蛋白)

- 中文名称: ACTR2 抗原(重组蛋白)
- 英文名称: ACTR2 Antigen (Recombinant Protein)
- 别名: ARP2
- 储存: 冷冻(-20℃)
- 相关类别: 抗原
- 概 述:

Fusion protein corresponding to C terminal 300 amino acids of human ACTR2

技术规格:

Full name:	actin related protein 2
Synonyms:	ARP2
Swissprot:	P61160
Gene Accession:	BC014546
Purity:	>85%, as determined by Coomassie blue stained SDS-PAG E
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The specific function of this gene has not yet been deter mined; however, the protein it encodes is known to be a major constituent of the ARP2/3 complex. This complex is located at the cell surface and is essential to cell shape a nd motility through lamellipodial actin assembly and protr usion. Two transcript variants encoding different isoforms



have been found for this gene.