

ADAM15 抗原(重组蛋白)

中文名称: ADAM15 抗原(重组蛋白)

英文名称: ADAM15 Antigen (Recombinant Protein)

别 名: ADAM metallopeptidase domain 15; MDC15

储 存: 冷冻(-20℃)

相关类别: 抗原

概 述:

Fusion protein corresponding to a region derived from 213-414 amino acids of human ADAM15

技术规格:

Full name:	ADAM metallopeptidase domain 15
Synonyms:	MDC15
Swissprot:	Q13444
Gene Accession:	BC014566
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a member of the ADAM (a disin tegrin and metalloproteinase) protein family. ADAM family members a re type I transmembrane glycoproteins known to be involved in cell a dhesion and proteolytic ectodomain processing of cytokines and adhe sion molecules. This protein contains multiple functional domains including a zinc-binding metalloprotease domain, a disintegrin-like domain, as well as a EGF-like domain. Through its disintegrin-like domain, this protein specifically interacts with the integrin beta chain, beta 3. It



also interacts with Src family protein-tyrosine kinases in a phosphoryl ation-dependent manner, suggesting that this protein may function in cell-cell adhesion as well as in cellular signaling. Multiple alternatively spliced transcript variants encoding distinct isoforms have been observed.