

ADAM20 抗原(重组蛋白)

中文名称: ADAM20 抗原(重组蛋白)

英文名称: ADAM20 Antigen (Recombinant Protein)

别 名: ADAM metallopeptidase domain 20

储 存: 冷冻(-20℃)

相关类别: 抗原

概 述:

Fusion protein corresponding to a region derived from 211-400 amino acids of human ADAM20

技术规格:

Full name:	ADAM metallopeptidase domain 20
Swissprot:	O43506
Gene Accession:	BC025378
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the ADAM (a disintegrin and metal loprotease domain) family. Members of this family are membrane-an chored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biological processes involving c ell-cell and cell-matrix interactions, including fertilization, muscle dev elopment, and neurogenesis. The expression of this gene is testis-sp ecific.

