

ADAMTSL1 抗原（重组蛋白）

中文名称：ADAMTSL1 抗原（重组蛋白）

英文名称：ADAMTSL1 Antigen (Recombinant Protein)

别名：ADAMTS-like 1; C9orf94; PUNCTIN; ADAMTSR1; ADAMTSL-1

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 116-329 amino acids of human ADAMTSL1

技术规格：

Full name:	ADAMTS-like 1
Synonyms:	C9orf94; PUNCTIN; ADAMTSR1; ADAMTSL-1
Swissprot:	Q8N6G6
Gene Accession:	BC100788
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a secreted protein and member of the ADAMTS (a disintegrin and metalloproteinase with thrombospondin motif) family. This protein lacks the metalloproteinase and disintegrin-like domains, which are typical of the ADAMTS family, but contains other ADAMTS domains, including the thrombospondin type 1 motif. This protein may have important functions in the extracellular matrix. Alternative

Alternative splicing results in multiple transcript variants encoding distinct proteins.