

ADIPOQ 抗原（重组蛋白）

中文名称：ADIPOQ 抗原（重组蛋白）

英文名称：ADIPOQ Antigen (Recombinant Protein)

别名：adiponectin, C1Q and collagen domain containing; ACDC; ADPN; APM1; APM-1; GBP28; ACRP30; ADIPQTL1

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 19-244 amino acids of human ADIPOQ

技术规格：

Full name:	adiponectin, C1Q and collagen domain containing
Synonyms:	ACDC; ADPN; APM1; APM-1; GBP28; ACRP30; ADIPQTL1
Swissprot:	Q15848
Gene Accession:	BC096308
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene is expressed in adipose tissue exclusively. It encodes a protein with similarity to collagens X and VIII and complement factor C1q. The encoded protein circulates in the plasma and is involved with metabolic and hormonal processes. Mutations in this gene are associated with adiponectin deficiency. Multiple alternatively spliced variants, encoding t

he same protein, have been identified.