

AGR3 抗原(重组蛋白)

- 中文名称: AGR3 抗原(重组蛋白)
- 英文名称: AGR3 Antigen (Recombinant Protein)
- 别名: Anterior gradient 3 homolog (Xenopus laevis); AG3; HAG3; hAG-3; BCMP11;

PDIA18

- 储存: 冷冻(-20℃)
- 相关类别: 抗原
- 概 述:

Fusion protein corresponding to a region derived from 22-166 amino acids of human AGR3

技术规格:

Full name:	Anterior gradient 3 homolog (Xenopus laevis)
Synonyms:	AG3; HAG3; hAG-3; BCMP11; PDIA18
Swissprot:	Q8TD06
Gene Accession:	BC058284
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Anterior gradient protein 3 homolog is a protein that in humans is encoded by the AGR3 gene. The deduced 166-amino acid protein h as an N-terminal signal sequence. Fluorescence-tagged BCMP11 loca lized in secretory or endosome-like organelles in transfected breast carcinoma cell lines. Immunohistochemical analysis showed that expr ession of BCMP11 was restricted in normal tissues, with significant s



taining of only the epithelial lining of colonic mucosa. However, BC
MP11 was detected in 43 of 58 (74%) breast cancer tissues tested,
but not in normal breast ductal epithelial cells. In cancerous epitheli
al cells, BCMP11 localized to the cytoplasm.