

兔抗 TNFSF12 多克隆抗体

中文名称：兔抗 TNFSF12 多克隆抗体

英文名称： Anti-TNFSF12 rabbit polyclonal antibody

别名： APO3L; DR3LG; TWEAK

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： TNFSF12

宿主： Rabbit

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:

The protein encoded by this gene is a cytokine that belongs to the tumor necrosis factor (TNF) ligand family. This protein is a ligand for the FN14/TWEAKR receptor. This cytokine has overlapping signaling functions with TNF, but displays a much wider tissue distribution. This cytokine, which exists in both membrane-bound and secreted forms, can induce apoptosis via multiple pathways of cell death in a cell type-specific manner. This cytokine is also found to promote proliferation and migration of endothelial cells, and thus acts as a regulator of

	f angiogenesis.
Applications:	WB, IHC, IF
Name of antibody:	TNFSF12
Immunogen:	Synthesized peptide derived from internal of human TNF12.
Full name:	tumor necrosis factor (ligand) superfamily, member 12
Synonyms :	APO3L; DR3LG; TWEAK
SwissProt:	O43508
IHC positive control:	Human brain tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	27 kDa
WB Positive control:	COLO cells, Huvec cells and 293 cells lysates
WB Recommended dilution:	500-3000
IF positive control:	Hela cells
IF Recommend dilution:	100-500



