

## 小鼠抗 MTFR2 单克隆抗体

中文名称：小鼠抗 MTFR2 单克隆抗体

英文名称：Anti-MTFR2 mouse monoclonal antibody

别名：mitochondrial fission regulator 2; DUFD1; FAM54A

相关类别：一抗

储存：冷冻（-20℃）

宿主：Mouse

抗原：MTFR2

反应种属：Human, Dog, Rat, Monkey

标记物：Unconjugate

克隆类型：mouse monoclonal

### 技术规格

<b>Background:</b>	May play a role in mitochondrial aerobic respiration essentially in the testis. Can also promote mitochondrial fission (By similarity).
<b>Applications:</b>	WB
<b>Name of antibody:</b>	MTFR2
<b>Immunogen:</b>	Fusion protein of human MTFR2
<b>Full name:</b>	mitochondrial fission regulator 2
<b>Synonyms:</b>	DUFD1; FAM54A
<b>SwissProt:</b>	Q6P444
<b>WB Predicted band size:</b>	43 KD
<b>WB Positive control:</b>	HepG2, Hela, SVT2, A549, COS7, Jurkat, MDCK, PC12 cell I



	ysates
<b>WB Recommended dilution:</b>	500-2000