

小鼠抗 ROMO1 单克隆抗体

中文名称: 小鼠抗 ROMO1 单克隆抗体

英文名称: Anti-ROMO1 mouse monoclonal antibody

别名: bA353C18.2; C20orf52; MTGM; MTGMP

储存: 冷冻 (-20℃) 避光

抗原: ROMO1

宿主: Mouse

反应种属: Human

标记物: Unconjugate

相关类别: 一抗

克隆类型: mouse monoclonal

技术规格

Background:	The protein encoded by this gene is a mitochondrial membrane protein that is responsible for increasing the level of reactive oxygen species (ROS) in cells. The protein also has antimicrobial activity against a variety of bacteria by inducing bacterial membrane breakage.
Applications:	WB, IHC, IF
Name of antibody:	ROMO1
Immunogen:	Fusion protein of human ROMO1
Full name:	reactive oxygen species modulator 1 (ROMO1)
Synonyms:	bA353C18.2; C20orf52; MTGM; MTGMP

SwissProt:	P60602
IHC positive control:	adenocarcinoma of human ovary tissue and adenocarcinoma of human endometrium tissue; carcinoma of human thyroid tissue and human breast tissue
IHC Recommend dilution:	30-150
WB Predicted band size:	8 kDa
WB Positive control:	HepG2, HeLa, SVT2, A549, COS7, Jurkat, MDCK cell lysates
WB Recommended dilution:	200-1000