

## 小鼠抗 NCOA7 单克隆抗体

中文名称： 小鼠抗 NCOA7 单克隆抗体

英文名称： Anti-NCOA7 mouse monoclonal antibody

别名： nuclear receptor coactivator 7; ESNA1; TLDC4; ERAP140; NCOA7-AS; Nbla00052; Nbla10993; dJ187J11.3

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Mouse

抗原： NCOA7

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	Enhances the transcriptional activities of several nuclear receptors. Involved in the coactivation of different nuclear receptors, such as ESR1, THRB, PPARG and RARA.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	NCOA7
<b>Immunogen:</b>	Fusion protein of human NCOA7
<b>Full name:</b>	nuclear receptor coactivator 7
<b>Synonyms:</b>	ESNA1; TLDC4; ERAP140; NCOA7-AS; Nbla00052; Nbla10993; dJ187J11.3
<b>SwissProt:</b>	Q8NI08



<b>WB Predicted band size:</b>	106 KD
<b>WB Positive control:</b>	K562, M7, Hela cell lysates
<b>WB Recommended dilution:</b>	500-2000