

## 小鼠抗 RUNX1T1 单克隆抗体

- 中文名称: 小鼠抗 RUNX1T1 单克隆抗体
- 英文名称: Anti-RUNX1T1 mouse monoclonal antibody
- 别 名: CDR; ETO; MTG8; AML1T1; ZMYND2; CBFA2T1; AML1-MTG8
- 储 存:冷冻(-20℃) 避光
- 抗 原: RUNX1T1
- 宿 主: Mouse
- 反应种属: Human
- 相关类别: 一抗
- 标记物: Unconjugate
- 克隆类型: mouse monoclonal

## 技术规格

Background:	This gene encodes a member of the myeloid translocation gene family which interact with DNA-bound transcription factors and r ecruit a range of corepressors to facilitate transcriptional repressi on. The t(8;21)(q22;q22) translocation is one of the most frequen t karyotypic abnormalities in acute myeloid leukemia. The translo cation produces a chimeric gene made up of the 5'-region of th e runt-related transcription factor 1 gene fused to the 3'-region of this gene. The chimeric protein is thought to associate with t he nuclear corepressor/histone deacetylase complex to block he matopoietic differentiation. Alternative splicing results in multiple
-------------	--



	transcript variants. [provided by RefSeq].
Applications:	WB, IF
Name of antibody:	RUNX1T1
Immunogen:	Fusion protein of human RUNX1T1
Full name:	runt-related transcription factor 1; translocated to, 1 (cyclin D-rel ated) (RUNX1T1), transcript variant 3
Synonyms:	CDR; ETO; MTG8; AML1T1; ZMYND2; CBFA2T1; AML1-MTG8
SwissProt:	Q06455
WB Predicted band siz e:	68 kDa
WB Positive control:	Hela, SVT2, A549, COS7, Jurkat, PC12, MCF-7 cell lysates
WB Recommended dilu tion:	500-2000