

## 小鼠抗 NEUROD4 单克隆抗体

中文名称： 小鼠抗 NEUROD4 单克隆抗体

英文名称： Anti-NEUROD4 mouse monoclonal antibody

别名： neuronal differentiation 4; ATH3; ATH-3; Atoh3; MATH3; MATH-3; bHLHa4

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： NEUROD4

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	Probably acts as a transcriptional activator. Mediates neuronal differentiation. Required for the regulation of amacrine cell fate specification in the retina (By similarity).
<b>Applications:</b>	WB
<b>Name of antibody:</b>	NEUROD4
<b>Immunogen:</b>	Fusion protein of human NEUROD4
<b>Full name:</b>	neuronal differentiation 4
<b>Synonyms:</b>	ATH3; ATH-3; Atoh3; MATH3; MATH-3; bHLHa4
<b>SwissProt:</b>	Q9HD90
<b>WB Predicted band size:</b>	37 KD



<b>WB Positive control:</b>	HepG2, Jurkat and MCF7 cell lysates
<b>WB Recommended dilution:</b>	200-1000