

小鼠抗 NFIC 单克隆抗体

中文名称： 小鼠抗 NFIC 单克隆抗体

英文名称： Anti-NFIC mouse monoclonal antibody

别名： nuclear factor I C; CTF; NFI; CTF5; NF-I

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： NFIC

反应种属： Human

标记物： Unconjugate

克隆类型： Mouse Monoclonal

技术规格

Background:	The protein encoded by this gene belongs to the CTF/NF-I family. These are dimeric DNA-binding proteins, and function as cellular transcription factors and as replication factors for adenovirus DNA replication. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.
Applications:	WB, IF
Name of antibody:	NFIC
Immunogen:	Fusion protein of human NFIC
Full name:	nuclear factor I C
Synonyms:	CTF; NFI; CTF5; NF-I

SwissProt:	P08651
WB Predicted band size:	56 kDa
WB Positive control:	Hela NE cell lysate
WB Recommended dilution:	200-1000
IF positive control:	Hela cells
IF Recommend dilution:	100-200



