

兔抗 TOLLIP 多克隆抗体

中文名称：兔抗 TOLLIP 多克隆抗体

英文名称：Anti-TOLLIP rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）避光

宿 主：Rabbit

抗 原：TOLLIP

反应种属：Human

标 记 物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Component of the signaling pathway of IL-1 and Toll-like receptors. Inhibits cell activation by microbial products. Recruits IRAK1 to the IL-1 receptor complex. Inhibits IRAK1 phosphorylation and kinase activity.
Applications:	WB
Name of antibody:	TOLLIP
Immunogen:	Synthesized peptide derived from Internal of human TOLLIP.
Full name:	toll interacting protein
Synonyms :	0
SwissProt:	Q9H0E2
WB Predicted band size:	30 kDa
WB Positive control:	Huvec, Jurkat and 293 cells lysates

WB Recommended dilution: 500-3000

