

小鼠抗 SERPINB8 单克隆抗体

中文名称: 小鼠抗 SERPINB8 单克隆抗体

英文名称: Anti-SERPINB8 mouse monoclonal antibody

别 名: serpin family B member 8; PI8; CAP2; HMSD; PI-8; PSS5; HSMD-v; C18orf53

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Mouse

抗 原: SERPINB8

反应种属: Human

标 记 物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	The protein encoded by this gene is a member of the ov-s erpin family of serine protease inhibitors. The encoded prot ein is produced by platelets and can bind to and inhibit the function of furin, a serine protease involved in platelet functions. In addition, this protein has been found to enhance the mechanical stability of cell-cell adhesion in the skin, and defects in this gene have been associated with an autosomal recessive form of extelliative inhibitors.
Applications:	osomal-recessive form of exfoliative ichthyosis. WB
Name of antibody:	SERPINB8



Immunogen:	Fusion protein of human SERPINB8
Full name:	serpin family B member 8
Synonyms:	PI8; CAP2; HMSD; PI-8; PSS5; HSMD-v; C18orf53
SwissProt:	P50452
WB Predicted band size:	43 KD
WB Positive control:	Jurkat, MCF7 cell lysates
WB Recommended dilution:	500-2000