

兔抗 TP53(Phospho-Ser20) 多克隆抗体

中文名称：兔抗 TP53(Phospho-Ser20) 多克隆抗体

英文名称： Anti-TP53(Phospho-Ser20) rabbit polyclonal antibody

别名： P53; BCC7; LFS1; TRP53

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： TP53(Phospho-Ser20)

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene encodes a tumor suppressor protein containing transcriptional activation, DNA binding, and oligomerization domains. The encoded protein responds to diverse cellular stresses to regulate expression of target genes, thereby inducing cell cycle arrest, apoptosis, senescence, DNA repair, or changes in metabolism. Mutations in this gene are associated with a variety of human cancers, including hereditary cancers such as Li-Fraumeni syndrome. Alternative splicing of this gene and the use of alternate promoters result in multiple transcript variants and isoforms.

Applications:	WB
Name of antibody:	TP53(Phospho-Ser20)
Immunogen:	Peptide sequence around phosphorylation site of Serine 20 (T-F-S(p)-D-L) derived from Human p53.
Full name:	tumor protein p53
Synonyms :	P53; BCC7; LFS1; TRP53
SwissProt:	P04637
WB Predicted band size:	44 kDa
WB Positive control:	Mouse brain tissue lysates
WB Recommended dilution:	500-1000

