

兔抗 TP53(Phospho-Thr387) 多克隆抗体

中文名称：兔抗 TP53(Phospho-Thr387) 多克隆抗体

英文名称： Anti-TP53(Phospho-Thr387) rabbit polyclonal antibody

别名： P53; BCC7; LFS1; TRP53

抗原： TP53(Phospho-Thr387)

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a tumor suppressor protein containing transcriptional activation, DNA binding, and oligomerization domains. The encoded protein responds to diverse cellular stresses to regulate expression of target genes, thereby inducing cell cycle arrest, apoptosis, senescence, DNA repair, or changes in metabolism. Mutations in this gene are associated with a variety of human cancers, including hereditary cancers such as Li-Fraumeni syndrome.
Applications:	IHC
Name of antibody:	TP53(Phospho-Thr387)

Immunogen:	Peptide sequence around phosphorylation site of threonine 387(F-K-T(p)-E-G) derived from Human p53.
Full name:	tumor protein p53
Synonyms :	P53; BCC7; LFS1; TRP53
SwissProt:	P04637
IHC positive control:	Human breast carcinoma tissue
IHC Recommend dilution:	50-100

