

小鼠抗 SIRPB1 单克隆抗体

中文名称: 小鼠抗 SIRPB1 单克隆抗体

英文名称: Anti-SIRPB1 mouse monoclonal antibody

别 名: signal regulatory protein beta 1; CD172b; SIRP-BETA-1

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Mouse

抗原: SIRPB1

反应种属: Human

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格



Name of antibody:	SIRPB1
Immunogen:	Fusion protein of human SIRPB1
Full name:	signal regulatory protein beta 1
Synonyms:	CD172b; SIRP-BETA-1
SwissProt:	O00241
WB Predicted band size :	43 KD
WB Positive control:	A549, COS7, MDCK cell lysates
WB Recommended dilut ion:	200-1000