

## 小鼠抗 SMC1A(N-term)单克隆抗体

中文名称: 小鼠抗 SMC1A(N-term)单克隆抗体

英文名称: Anti-SMC1A(N-term) mouse monoclonal antibody

别 名: structural maintenance of chromosomes 1A; SMC1; SMCB; CDLS2; SB1.8; SMC1L1; DXS423E; SMC1alpha

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Mouse

抗 原: SMC1A(N-term)

反应种属: Human

标 记 物: Unconjugate

克隆类型: Mouse Monoclonal

## 技术规格

**Background:** 

Proper cohesion of sister chromatids is a prerequisite for the correc t segregation of chromosomes during cell division. The cohesin mul tiprotein complex is required for sister chromatid cohesion. This complex is composed partly of two structural maintenance of chromo somes (SMC) proteins, SMC3 and either SMC1B or the protein encoded by this gene. Most of the cohesin complexes dissociate from the chromosomes before mitosis, although those complexes at the kinetochore remain. Therefore, the encoded protein is thought to be an important part of functional kinetochores. In addition, this protein interacts with BRCA1 and is phosphorylated by ATM, indicating a



	potential role for this protein in DNA repair. This gene, which belon gs to the SMC gene family, is located in an area of the X-chromos ome that escapes X inactivation. Mutations in this gene result in Co rnelia de Lange syndrome. Alternative splicing results in multiple tra nscript variants encoding different isoforms.
Applications:	WB, IHC, FC
Name of antibody:	SMC1A(N-term)
Immunogen:	Fusion protein of human SMC1A
Full name:	structural maintenance of chromosomes 1A
Synonyms:	SMC1; SMCB; CDLS2; SB1.8; SMC1L1; DXS423E; SMC1alpha
SwissProt:	Q14683
IHC positive control:	Human colon
IHC Recommend dil ution:	500-1000
WB Predicted band size:	143 kDa
WB Positive control:	HelaNE, Jurkat and A431 cell lysates
WB Recommended dilution:	1000-5000
FC Recommend dilu tion:	50-100









