

## 兔抗 TRA2A 多克隆抗体

中文名称：兔抗 TRA2A 多克隆抗体

英文名称： Anti-TRA2A rabbit polyclonal antibody

别名： transformer 2 alpha homolog (Drosophila)

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： TRA2A

宿主： Rabbit

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is a member of the transformer 2 homolog family and encodes a protein with several RRM (RNA recognition motif) domains. This phosphorylated nuclear protein binds to specific RNA sequences and plays a role in the regulation of pre-mRNA splicing. Alternatively spliced transcript variants have been found for this gene.
<b>Applications:</b>	WB, IHC, IF
<b>Name of antibody:</b>	TRA2A
<b>Immunogen:</b>	Synthesized peptide derived from internal of human TRA-2 $\alpha$ .
<b>Full name:</b>	transformer 2 alpha homolog

<b>Synonyms:</b>	AWMS1; HSU53209
<b>SwissProt:</b>	Q13595
<b>IHC positive control:</b>	Human brain tissue
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	33 kDa
<b>WB Positive control:</b>	Huvec cells lysate
<b>WB Recommended dilution:</b>	500-3000
<b>IF positive control:</b>	Hela cells
<b>IF Recommend dilution:</b>	100-500



