

小鼠抗 TUBB2B 单克隆抗体

中文名称: 小鼠抗 TUBB2B 单克隆抗体

英文名称: Anti-TUBB2B mouse monoclonal antibody

别名: bA506K6.1; PMGYSA

抗原: TUBB2B

储存: 冷冻 (-20℃) 避光

宿主: Mouse

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	The protein encoded by this gene is a beta isoform of tubulin, which binds GTP and is a major component of microtubules. This gene is highly similar to TUBB2A and TUBB2C. Defects in this gene are a cause of asymmetric polymicrogyria.
Applications:	WB, IHC, IF
Name of antibody:	TUBB2B
Immunogen:	Fusion protein of human TUBB2B
Full name:	tubulin, beta 2B class IIb (TUBB2B)
Synonyms:	bA506K6.1; PMGYSA
SwissProt:	Q9BVA1
IHC positive control:	adenocarcinoma of human ovary tissue and human kidney tissue

IHC Recommend dilution:	25-100
WB Predicted band size:	50 kDa
WB Positive control:	HepG2, Hela, SVT2, COS7, MDCK, PC12 cell lysates
WB Recommended dilution:	500-2000