

## 小鼠抗 UBE2S 单克隆抗体

中文名称: 小鼠抗 UBE2S 单克隆抗体

英文名称: Anti-UBE2S mouse monoclonal antibody

别名: E2-EPF; E2EPF; EPF5

抗原: UBE2S

储存: 冷冻 (-20℃) 避光

宿主: Mouse

相关类别: 一抗

反应种属: Human , Monkey , Mouse

标记物: Unconjugate

克隆类型: mouse monoclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the ubiquitin-conjugating enzyme family. The encoded protein is able to form a thiol ester linkage with ubiquitin in a ubiquitin activating enzyme-dependent manner, a characteristic property of ubiquitin carrier proteins. [provided by RefSeq, Jul 2008].
<b>Applications:</b>	WB, IF, FC
<b>Name of antibody:</b>	UBE2S
<b>Immunogen:</b>	Fusion protein of human UBE2S
<b>Full name:</b>	ubiquitin-conjugating enzyme E2S (UBE2S)
<b>Synonyms:</b>	E2-EPF; E2EPF; EPF5
<b>SwissProt:</b>	Q16763
<b>WB Predicted band size:</b>	24 kDa
<b>WB Positive control:</b>	HepG2, HeLa, SVT2, A549, COS7, Jurkat, MCF-7 cell lysates



<b>WB Recommended dilution:</b>	200-1000
---------------------------------	----------