

小鼠抗 VPS28 单克隆抗体

中文名称: 小鼠抗 VPS28 单克隆抗体

英文名称: Anti-VPS28 mouse monoclonal antibody

别名: MGC60323

抗原: VPS28

储存: 冷冻 (-20°C) 避光

宿主: Mouse

反应种属: Human , Dog , Rat , Monkey , Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	This gene encodes a protein subunit of the ESCRT-I complex (endosomal complexes required for transport), which functions in the transport and sorting of proteins into subcellular vesicles. This complex can also be hijacked to facilitate the budding of enveloped viruses from the cell membrane. Alternative splicing results in multiple transcript variants encoding different isoforms.
Applications:	WB, IHC, IF
Name of antibody:	VPS28
Immunogen:	Fusion protein of human VPS28
Full name:	VPS28 subunit of ESCRT-I
SwissProt:	Q9UK41
IHC positive control:	human pancreas tissue

IHC Recommend dilution:	30-150
WB Predicted band size:	25 kDa
WB Positive control:	HepG2, Hela, SVT2, A549, COS7, Jurkat, MDCK, PC12, MCF-7 cell lysates
WB Recommended dilution:	200-1000
IF positive control:	COS7 cells transiently transfected by pCMV6-ENTRY
IF Recommend dilution:	100-200