

## 小鼠抗 VPS50 单克隆抗体

中文名称：小鼠抗 VPS50 单克隆抗体

英文名称：Anti-VPS50 mouse monoclonal antibody

别名：VPS50, EARP/GARPII complex subunit; VPS54L; CCDC132

相关类别：一抗

储存：冷冻（-20℃）

宿主：Mouse

抗原：VPS50

反应种属：Human

标记物：Unconjugate

克隆类型：mouse monoclonal

### 技术规格

<b>Background:</b>	Acts as component of the EARP complex that is involved in endocytic recycling. The EARP complex associates with Rab4-positive endosomes and promotes recycling of internalized transferrin receptor (TFRC) to the plasma membrane. Within the EARP complex, required to tether the complex to recycling endosomes. Not involved in retrograde transport from early and late endosomes to the trans-Golgi network (TGN).
<b>Applications:</b>	WB
<b>Name of antibody:</b>	VPS50
<b>Immunogen:</b>	Fusion protein of human VPS50
<b>Full name:</b>	VPS50, EARP/GARPII complex subunit
<b>Synonyms:</b>	VPS54L; CCDC132

<b>SwissProt:</b>	Q96JG6
<b>WB Predicted band size:</b>	111 KD
<b>WB Positive control:</b>	Hela cell lysate
<b>WB Recommended dilution:</b>	200-1000