

小鼠抗 WNT3A 单克隆抗体

中文名称: 小鼠抗 WNT3A 单克隆抗体

英文名称: Anti-WNT3A mouse monoclonal antibody

别 名: Wnt family member 3A

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Mouse

抗 原: WNT3A

反应种属: Human

标 记 物: Unconjugate

克隆类型: Mouse Monoclonal

技术规格

Background:	The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family. It encodes a protein which shows 96% amino acid identity to mouse Wnt3A protein, and 84% to human WNT3 protein, another WNT gene product. This gene is clustered with WNT14 gene, another family member, in chromosome 1q42 region.
Applications:	WB
Name of antibody:	WNT3A
Immunogen:	Fusion protein of human WNT3A



Full name:	Wnt family member 3A
SwissProt:	P56704
WB Predicted band size:	39 kDa
WB Positive control:	MCF-7, ZR75-1, MDA-MB-468 and T47D cell lysates
WB Recommended dilution:	500-2000

