

## RPS11 抗原(重组蛋白)

中文名称: RPS11 抗原(重组蛋白)

英文名称: RPS11 Antigen (Recombinant Protein)

别 名: ribosomal protein S11; S11

储 存: 冷冻(-20℃)

相关类别: 抗原

## 概述

Fusion protein corresponding to a region derived from 2-158 amino acids of human RPS11

## 技术规格

Full name:	ribosomal protein S11
Synonyms:	S11
Swissprot:	P62280
Gene Accession:	BC007283
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 struct urally distinct proteins. This gene encodes a member of the S17P f amily of ribosomal proteins that is a component of the 40S subunit. This gene is co-transcribed with the small nucleolar RNA gene U 35B, which is located in the third intron. As is typical for genes en



coding ribosomal proteins, there are multiple processed pseudogen	
es of this gene dispersed throughout the genome.	