

THAP4 抗原（重组蛋白）

中文名称： THAP4 抗原（重组蛋白）

英文名称： THAP4 Antigen (Recombinant Protein)

别名： THAP domain containing 4; PP238; CGI-36

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to C terminal 252 amino acids of human THAP4

技术规格

Full name:	THAP domain containing 4
Synonyms:	PP238; CGI-36
Swissprot:	Q8WY91
Gene Accession:	BC000247
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	THAP domain-containing protein 4 is a protein that in humans is encoded by the THAP4 gene.