

TLE1 抗原（重组蛋白）

中文名称： TLE1 抗原（重组蛋白）

英文名称： TLE1 Antigen (Recombinant Protein)

别名： transducin like enhancer of split 1; ESG; ESG1; GRG1

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 571-770 amino acids of human TLE1

技术规格

Full name:	transducin like enhancer of split 1
Synonyms:	ESG; ESG1; GRG1
Swissprot:	Q04724
Gene Accession:	BC015747
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Transcriptional corepressor that binds to a number of transcription factors. Inhibits NF-kappa-B-regulated gene expression. Inhibits the transcriptional activation mediated by FOXA2, and by CTNNB1 and TCF family members in Wnt signaling. The effects of full-length TLE family members may be modulated by association with dominant-negative AES. Unusual function as coactivator for ESRRG.