

## TMEFF1 抗原（重组蛋白）

中文名称： TMEFF1 抗原（重组蛋白）

英文名称： TMEFF1 Antigen (Recombinant Protein)

别名： transmembrane protein with EGF like and two follistatin like domains 1; TR-1; H7365; C9orf2; CT120.1

储存： 冷冻（-20℃）

相关类别： 抗原

### 概述

Fusion protein corresponding to a region derived from 51-250 amino acids of human TMEFF1

### 技术规格

<b>Full name:</b>	transmembrane protein with EGF like and two follistatin like domains 1
<b>Synonyms:</b>	TR-1; H7365; C9orf2; CT120.1
<b>Swissprot:</b>	Q8IYR6
<b>Gene Accession:</b>	BC035056
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	May inhibit NODAL and BMP signaling during neural patterning (By similarity). May be a tumor suppressor in brain cancers.