

TMPRSS11A 抗原（重组蛋白）

中文名称： TMPRSS11A 抗原（重组蛋白）

英文名称： TMPRSS11A Antigen (Recombinant Protein)

别名： transmembrane protease, serine 11A; ECRG1

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 203-353 amino acids of human TMPRSS11A

技术规格

Full name:	transmembrane protease, serine 11A
Synonyms:	ECRG1
Swissprot:	Q6ZMR5
Gene Accession:	BC111796
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Probable serine protease which may play a role in cellular senescence. Overexpression inhibits cell growth and induce G1 cell cycle arrest.