

VPS33A 抗原(重组蛋白)

- 中文名称: VPS33A 抗原(重组蛋白)
- 英文名称: VPS33A Antigen (Recombinant Protein)
- 别名: VPS33A core subunit of CORVET and HOPS complexes; MPSPS
- 储存: 冷冻(-20℃)
- 相关类别: 抗原

概述

Fusion protein corresponding to C terminal 200 amino acids of human VPS33A

技术规格

Full name:	VPS33A core subunit of CORVET and HOPS complexes
Synonyms:	MPSPS
Swissprot:	Q96AX1
Gene Accession:	BC016617
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a tethering protein and a core subunit of the homotypic fusion and protein sorting (HOPS) complex. Th e HOPS complex and a second endosomal tethering complex called the class C core vacuole/endosome tethering (CORVET) complex, perform diverse functions in endocytosis including m embrane tethering, RabGTPase interaction, activation and proo freading of synaptic-soluble N-ethylmaleimide-sensitive factor



attachment receptor (SNARE) assembly to drive membrane fus ion, and endosome-to-cytoskeleton attachment. The HOPS co mplex controls endosome maturation as well as endosome tra ffic to the lysosome. This complex is essential for vacuolar fus ion and is required for adaptor protein complex 3-dependent transport from the golgi to the vacuole. The encoded protein belongs to the Sec1/Munc18 (SM) family of SNARE-mediated membrane fusion regulators. Naturally occurring mutations in this gene are associated with a novel mucopolysaccharidosis-li ke disease.