

## 兔抗 ENPP2 多克隆抗体

中文名称:	兔抗 ENPP2 多克隆抗体
英文名称:	Anti-ENPP2 rabbit polyclonal antibody
别名:	ATX; NPP2; ATX-X; PDNP2; LysoPLD; AUTOTAXIN; PD-IALPHA
相关类别:	一抗
储存:	冷冻 (-20℃)
宿主:	Rabbit
抗原:	ENPP2
反应种属:	Human, Mouse, Rat
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene functions as both a phosphodiesterase, which cleaves phosphodiester bonds at the 5' end of oligonucleotides, and a phospholipase, which catalyzes production of lysophosphatidic acid (LPA) in extracellular fluids. LPA evokes growth factor-like responses including stimulation of cell proliferation and chemotaxis. This gene product stimulates the motility of tumor cells and has angiogenic properties, and its expression is upregulated in several kinds of carcinomas. The gene product is secreted and further processed to make the biologically active form. Several alternatively spliced transcript variants encoding different isoforms have been identified.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ENPP2
<b>Immunogen:</b>	Synthetic peptide of human ENPP2

<b>Full name:</b>	ectonucleotide pyrophosphatase/phosphodiesterase 2
<b>Synonyms:</b>	ATX; NPP2; ATX-X; PDNP2; LysoPLD; AUTOTAXIN; PD-IALP HA
<b>SwissProt:</b>	Q13822
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer and human brain
<b>IHC Recommend dilution:</b>	30-150



