

小鼠抗 NFKBIA(N-term)单克隆抗体

- 中文名称: 小鼠抗 NFKBIA(N-term)单克隆抗体
- 英文名称: Anti-NFKBIA(N-term) mouse monoclonal antibody
- 别名: NFKB inhibitor alpha; IKBA; MAD-3; NFKBI
- 相关类别: 一抗
- 储存: 冷冻 (-20℃)
- 宿主: Mouse
- 抗原: NFKBIA(N-term)
- 反应种属: Human, Mouse
- 标记物: Unconjugate
- 克隆类型: Mouse Monoclonal

技术规格

Background:	This gene encodes a member of the NF-kappa-B inhibitor family, which contain multiple ankrin repeat domains. The encoded protein interacts with REL dimers to inhibit NF-kappa-B/REL complexes which are involved in inflammatory responses. The encoded protein moves between the cytoplasm and the nucleus via a nuclear localization signal and CRM1-mediated nuclear export. Mutations in this gene have been found in ectodermal dysplasia anhidrotic with T-cell immunodeficiency autosomal dominant disease.
Applications:	WB
Name of antibody:	NFKBIA(N-term)

Immunogen:	Fusion protein of human NFKBIA
Full name:	NFKB inhibitor alpha
Synonyms:	IKBA; MAD-3; NFKBI
SwissProt:	P25963
WB Predicted band size:	36 kDa
WB Positive control:	Jurkat and Hela cell lysates
WB Recommended dilution:	1000-5000

