

## 兔抗 TOE1 多克隆抗体

中文名称：兔抗 TOE1 多克隆抗体

英文名称：Anti-TOE1 rabbit polyclonal antibody

别名：target of EGR1, exonuclease; PCH7; hCaf1z

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：TOE1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Inhibits cell growth rate and cell cycle. Induces CDKN1A expression as well as TGF-beta expression. Mediates the inhibitory growth effect of EGR1. Involved in the maturation of snRNAs and snRNA 3'-tail processing (PubMed:28092684).
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	TOE1
<b>Immunogen:</b>	Fusion protein of human TOE1
<b>Full name:</b>	target of EGR1, exonuclease
<b>Synonyms:</b>	PCH7; hCaf1z
<b>SwissProt:</b>	Q96GM8
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil

<b>IHC Recommend dilution:</b>	100-200
<b>WB Predicted band size:</b>	57 kDa
<b>WB Positive control:</b>	K562 cell lysate
<b>WB Recommended dilution:</b>	500-2000



