

兔抗 TPMT 多克隆抗体

中文名称: 兔抗 TPMT 多克隆抗体

英文名称: Anti-TPMT rabbit polyclonal antibody

别 名: thiopurine S-methyltransferase; TPMTD

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Rabbit

抗 原: TPMT

反应种属: Human

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes the enzyme that metabolizes thiopur ine drugs via S-adenosyl-L-methionine as the S-methyl donor and S-adenosyl-L-homocysteine as a byproduct. Thiopurine drugs such as 6-mercaptopurine are used as chemotherapeutic agents. Genetic polymorphisms that a ffect this enzymatic activity are correlated with variation s in sensitivity and toxicity to such drugs within individ uals, causing thiopurine S-methyltransferase deficiency. Related pseudogenes have been identified on chromoso mes 3, 18 and X.
Applications:	ELISA, IHC
Name of antibody:	TPMT



Immunogen:	Fusion protein of human TPMT
Full name:	thiopurine S-methyltransferase
Synonyms:	TPMTD
SwissProt:	P51580
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	40-200

