

兔抗 TRADD 多克隆抗体

中文名称：兔抗 TRADD 多克隆抗体

英文名称： Anti-TRADD rabbit polyclonal antibody

别名： Hs.89862

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： TRADD

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a death domain containing adaptor molecule that interacts with TNFRSF1A/TNFR1 and mediates programmed cell death signaling and NF-kappaB activation. This protein binds adaptor protein TRAF2, reduces the recruitment of inhibitor-of-apoptosis proteins (IAPs) by TRAF2, and thus suppresses TRAF2 mediated apoptosis. This protein can also interact with receptor TNFRSF6/FAS and adaptor protein FADD/MORT1, and is involved in the Fas-induced cell death pathway.
Applications:	ELISA, IHC
Name of antibody:	TRADD
Immunogen:	Fusion protein of human TRADD

Full name:	TNFRSF1A associated via death domain
Synonyms:	Hs.89862
SwissProt:	Q15628
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and Human thyroid cancer
IHC Recommend dilution:	100-300



