

## 兔抗 TRAPPC4 多克隆抗体

中文名称：兔抗 TRAPPC4 多克隆抗体

英文名称：Anti-TRAPPC4 rabbit polyclonal antibody

别名：SBDN; TRS23; PTD009; CGI-104; HSPC172; SYNBINDIN

相关类别：一抗

抗原：TRAPPC4

储存：冷冻 (-20℃)

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	May play a role in vesicular transport from endoplasmic reticulum to Golgi.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	TRAPPC4
<b>Immunogen:</b>	Fusion protein of human TRAPPC4
<b>Full name:</b>	trafficking protein particle complex 4
<b>Synonyms:</b>	SBDN; TRS23; PTD009; CGI-104; HSPC172; SYNBINDIN
<b>SwissProt:</b>	Q9Y296
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer and Human thyroid cancer

IHC Recommend dilution:

100-300



