

兔抗 ANGPTL6 多克隆抗体

中文名称：兔抗 ANGPTL6 多克隆抗体

英文名称：Anti-ANGPTL6 rabbit polyclonal antibody

别名：AGF; ARP5

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ANGPTL6

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May play a role in the wound healing process. May promote epidermal proliferation, remodeling and regeneration. May promote the chemotactic activity of endothelial cells and induce neovascularization. May counteract high-fat diet-induced obesity and related insulin resistance through increased energy expenditure.
Applications:	ELISA, IHC
Name of antibody:	ANGPTL6
Immunogen:	Synthetic peptide of human ANGPTL6
Full name:	angiotensin like 6
Synonyms:	AGF; ARP5
SwissProt:	Q8NI99

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	25-100

