

兔抗 ANXA4 多克隆抗体

- 中文名称: 兔抗 ANXA4 多克隆抗体
- 英文名称: Anti-ANXA4 rabbit polyclonal antibody
- 别名: ANX4; PIG28; ZAP36
- 抗原: ANXA4
- 储存: 冷冻(-20℃)
- 宿 主: Rabbit
- 相关类别: 一抗
- 反应种属: Human, Mouse, Rat
- 标记物: Unconjugate
- 克隆类型: rabbit polyclonal

技术规格

| Background: | Annexin IV (ANX4) belongs to the annexin family of calciu m-dependent phospholipid binding proteins. Although thei r functions are still not clearly defined, several members o f the annexin family have been implicated in membrane-r elated events along exocytotic and endocytotic pathways. ANX4 has 45 to 59% identity with other members of its f amily and shares a similar size and exon-intron organizati on. Isolated from human placenta, ANX4 encodes a protei n that has possible interactions with ATP, and has in vitro anticoagulant activity and also inhibits phospholipase A2 a ctivity. ANX4 is almost exclusively expressed in epithelial c |
|-------------|---|
| | ctivity. ANX4 is almost exclusively expressed in epithelial c ells. |



| Applications: | ELISA, WB, IHC |
|-----------------------------|---|
| Name of antibody: | ANXA4 |
| Immunogen: | Fusion protein of human ANXA4 |
| Full name: | annexin A4 |
| Synonyms : | ANX4; PIG28; ZAP36 |
| SwissProt: | P09525 |
| ELISA Recommended dilution: | 1000-2000 |
| IHC positive control: | Human thyroid cancer |
| IHC Recommend dilution: | 15-50 |
| WB Recommended dilution: | 200-1000 |
| WB Positive control: | Raw264.7 cells and human fetal liver tissue |
| WB Predicted band size: | 36 kDa |





