

## 兔抗 APCDD1 多克隆抗体

中文名称：兔抗 APCDD1 多克隆抗体

英文名称：Anti-APCDD1 rabbit polyclonal antibody

别名：HHS; HTS; B7323; HYPT1; DRAPC1; FP7019

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：APCDD1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This locus encodes an inhibitor of the Wnt signaling pathway. Mutations at this locus have been associated with hereditary hypotrichosis simplex. Increased expression of this gene may also be associated with colorectal carcinogenesis.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	APCDD1
<b>Immunogen:</b>	Synthetic peptide of human APCDD1
<b>Full name:</b>	APC down-regulated 1
<b>Synonyms:</b>	HHS; HTS; B7323; HYPT1; DRAPC1; FP7019
<b>SwissProt:</b>	Q8J025
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human colorectal cancer
<b>IHC Recommend dilution:</b>	50-100

